

AIRBRUSH MAINTENANCE PROCEDURE

There is required maintenance to effectively maintain a well functioning airbrush applicator. Please follow the steps and images below after every use. The airbrush applicator must be cleaned **after every treatment use**. Failing to do so may jeopardize the warranty for this item.

The warranty for Oxygen Infusion Airbrush is 90-days. If the airbrush becomes compromised please follow the troubleshooting guidelines in order to resolve any concerns. If once you have followed those important steps and the airbrush is still impaired please call Skin for Life Service Dept. (866) 312.7546 for assistance.



STEP 1

- Take lid off from airbrush cup
- Pour O2 Cleaner into cup. Fill approximately 1/4 full
- Close cup with the lid



STEP 2

- Place a 4 x 4 gauze in one hand.
- Make sure tip of airbrush is touching the gauze. The liquid will build up and dislodge any solids
- Use index finger to push down and pull back the lever/trigger

STEP 2 continued

- Spray until the O2 cleaner is empty. (30 second process)

Step 3

- ADD TAP WATER - pour into the cup. Fill approximately 1/4 full
- Repeat the previous steps.

