Ultra Scrubber & Ultra Sonic Wave System

\$895.00

(Taxes & Shipping not included)

*Ask us about our monthly specials, promotion codes, or volume discounts.



Skin for ife sales@skinforlife.com


*prices may vary

www.skinforlife.com

www.skinforlifepurescience.com







Ultra Scrubber: Ultra Sonic Wave Esthetic Equipment

Dimensions: 11"L x 13.75"W x 4.75"H

Machine Weight: approximately 14 lbs.

Ultra Scrubber is fantastic for lifting dead skin cells without the friction or abrasion that may cause sensitivities to various skin conditions. Perfect for sun damaged skin as well as various pigmentation that are associated with sun.

Use the Ultra Scrubber spatula - blade for exfoliating, minimizing pigment irregularities, fine lines, and wrinkles. Ultra scrubber is also excellent for oily and troubled skin conditions, you may perform ultra-scrubber on mild, moderate, to severe oily and blemished skins.

The Buzz on Ultra Sonic Waves -

Using ultra sonic waves in the beauty industry, brought a heightened awareness when working with professional skin care formulas. The cavitation responses gave way to decreasing the thickness of stratum corneum, where professional treatment formulas have a greater chance of interacting with cell membranes.

Professionals are able to interact rapidly, and proficiently by delivering vital nutrients to skin. This is a great help when making skin change.

Ultra sonic waves through cavitation helps to lift, hydrate, and protect the skin's surface. Ultra Scrubber lifts by exfoliating the surface of the skin and ultra-sonic wave attachments enhance product penetration to support skin health.

Ultra Sonic waves are designed to increase the depths of product penetration. It allows skin professionals to indirectly affect deeper layers to encourage a healthier skin function. For mature skins, i.e. wrinkles, fine lines, and scars, sonic waves build volume - fullness for a beautiful youthful complexion.

System includes 3 handpieces:

Ultra Scrubber Hand Piece – spatula - blade for exfoliating, minimizing pigment irregularities, fine lines, and wrinkles.

Flat Ultra Sonic Hand Piece - 1 MHz hand piece (flat piece for the entire face)

The flat hand piece is for flat areas of the face where you may perform more lifting and firming movements to increase skin volume and fullness.

Cone Shape Ultra Sonic Hand Piece - 2 MHz hand piece (curved areas lip line, crow's feet)

This hand piece is a specialty hand piece that is used for the surrounding eye area or for more targeted areas of the face, i.e. around lip line because of fine lines or wrinkles.

Skin for Life Pure Science skin care included: (For Professional Use Only)

Enzyme Exfoliant + papaya complex with LSS™ delivery - 4oz.

Conductivity Gel - Hypoallergenic, Non-greasy, non-irritating, odorless and non-abrasive - 8oz.

Visit www.skinforlife.com for additional details with ingredients, protocols, and treatment suggestions.

Ultra Sonic Waves are excellent in conjunction with microdermabrasion, cosmetic acid exfoliation, microcurrent, oxygen infusion systems, and many more combinations for your clienteles' long term results.

*Prices subject to change without notice. Packaging subject to change.